

FACILIPRO™

Clean Dose Presoak Tablet



Enzymatic formula breaks down protein and starch films



Effective and efficient

- Presoaking flatware results in consistently clean and spotless results with every wash cycle
- Saves time and labor by helping to reduce need for rewash



User-friendly design

- Safer product handling
- Requires no PPE



Easy on the environment

- Concentrated solid tablets contribute to reduced water, energy, and plastic waste

The tablet's solid chemistry is formulated to save on packaging.

99.9%


reduction in plastic waste¹

The **Clean Dose Presoak Tablet** is part of the **Clean Dose Program**: an effective, easy-to-use tablet cleaning program available in a variety of applications to cover all your everyday cleaning needs.

1. Reduction in plastic waste is an approximation assuming use of one case of product per month over 12 months, compared with the equivalent use of Ecolab liquid products, which is enabled by the smaller, lighter package size for solid tablets.

ECOLAB®

A complete, innovative program that delivers a professional-grade clean for your customers and a clear, uncomplicated routine for your staff.



Presoak Tablet Directions for Use

- 1

Tear open a packet and drop the tablet into a presoak tub.

Wash hands after handling product.
- 2

Fill up the tub with warm water (110° – 120° F) and allow the tablet to dissolve.
- 3

Place dirty flatware into the solution and allow to soak.

Empower your employees with easy-to-use, color-coded tablets.

PRODUCT NAME	SKU	DILUTION RATIO	PACK SIZE
Clean Dose Presoak Tablet	6102799	1 tablet / tub	48 - 25 g



We thought small so you can think big.

Ready to learn more about our full tablet cleaning program?
Visit www.gofacilipro.com or scan the QR code.

